

● TOYOBO ENZYMES ●

(Industrial Grade)

IMMOBILIZED LIPASE

from Pseudomonas sp.

PREPARATION and SPECIFICATION

The enzyme is a preparation of a *Pseudomonas* sp. lipase (TOYOBO lipoprotein lipase Grade III LPL-311) immobilized on Hyflo Super-Cel.

Appearance	: Light brown powder (immobilized on Hyflo Super-Cel)
Activity	: 0.5U/mg-solid or more
Stabilizers	: Sugars

APPLICATIONS^{3,11)}

This enzyme is useful for enzymatic ester synthesis, transesterification, acidolysis and alcoholysis in organic solvents or solvent free substrate.

REFERENCES

- 1) *Biosci. Biotech. Biochem.*, 56 (7), 1118–1123 (1992)
- 2) *Agric. Biol. Chem.*, 53 (7), 1879–1884 (1989)
- 3) *J. Bull. Inst. Chem. Res., Kyoto Univ.* 70 (3), 333–337 (1992)
- 4) *Tetrahedron Letters*, 31, No.25, 3603–3604 (1990)
- 5) *Chemistry Letters*, 741–744 (1993)
- 6) *Tetrahedron Letters*, 35 No.43, 7997–8000 (1994)
- 7) *Tetrahedron: Asymmetry* 6, No.7, 1617–1620 (1995)
- 8) *Tetrahedron* 51 No.34, 9339–9352 (1995)
- 9) *Biosci. Biotech. Biochem.*, 59 (11), 2178–2180 (1995)
- 10) *Macromolecules* 28, No.21, 7046–7050 (1995)
- 11) *Journal of Fermentation and Bioengineering* 81 No.3, 206–211 (1996)

ASSAY

Principle:



Unit definition:

One unit causes the formation of the one micromole of ethyloctanoate (EO) per minute under the conditions described below.

Procedure

1. Pipette 10.0ml of dry diisopropyl ether containing 50mM 2,2,2-trifluoroethyl octanoate (TFEO) and 0.2M ethanol into a test tube with a magnetic stirrer bar and a stopper. Equilibrate at 25.5°C for about 10min with stirring.
2. Add about 20mg of immobilized lipase powder and stir constantly.
3. Take an aliquot of the assay mixture at a specific time interval and analyse with capillary gas chromatography (DB-5, J&W Scientific, 160°C). The retention time of TFEO and the product EO are 2.14 and 2.54 min, respectively.
4. A plot natural logarithm of the reaction conversion (C) vs reaction time gives a straight line $[-\ln(1-C) = k_0 t]$, and the initial rate (V_0) of the reaction is calculated from the slope of the line.

Calculation

Activity can be calculated by using the following formula :

$$\text{Initial rate (V}_0) = k_0 \times 50$$

$$\text{Weight activity (U/mg)} = \frac{V_0 \times 10}{W_s}$$

k_0 : Slope of the line $[-\ln(1-C) = k_0 t]$


50 : Concentration of TFEO (mM)

W_s : Sample weight (mg)


10 : Total volume (ml)

Illustrations of Stereoselective Acylation


Scheme1. Stereoselective Acylation of Binaphthol ²⁾



Scheme2. Stereoselective Acylation of Mandelonitrile ³⁾




Scheme3. Stereoselective Acylation of Tricarbonyl (o-trimethylsilylbenzyl alcohol) chromium ⁴⁾



活性測定法 (Japanese)

1.原理



2.定義

下記条件下で求めた、反応初期速度より算出した1分間に1マイクロモルのethyloctanoate(EO)を生成する酵素活性を1Uとする。

3.手順

- ①50mM 2,2,2-trifluoroethyl octanoate (TFEO)及び、0.2Mエタノールを含む乾燥diisopropyl ether [(i-Pr)₂O] 10mlをスターラーバーの入った共栓付試験管に採り、攪拌しながら25.5°Cで10分間予備加温する。
- ②固定化酵素 約20mgを上記反応液に添加し、1~20分間に数回サンプリングし、キャピラリーガスクロマトグラフィー (DB-5, J&W Scientific製, 160°C)にて分析する。TFEO及びEOのretention timeは 約2.14及び2.54分である。
- ③EOへの反応の転換率(C)より、反応時間t(分)に対して $-\ln(1-C)$ をプロットし、その傾き(k_0)より下記計算式に従って初期速度(V_0)及び比活性(U/mg)を求める。

4.計算式

$$V_0 = k_0 \times 50$$

$$U/mg = \frac{V_0 \times 10}{W_s}$$

- 50 : TFEO濃度 (mM)
 K_0 : 速度定数
 V_0 : 初期速度
 U/mg : 比活性
 W_s : サンプル重量 (mg)